

## Introduction

Invented by Binnig and Rohrer at IBM in 1981 (Nobel Prize in 1986).

Binnig also invented the Atomic Force Microscope (AFM) at Stanford University in 1986.

First instrument to give 3D

images of solid surface with atomic resolution.

Formation of nanostructures by localised heating or inducing chemical reactions.

Magnification 10

to 10  
9  
in x, y, and z directions  
with high resolution

## Theory and Principle

STM employs principle of electron tunneling.

The quantum mechanical equation assign a non zero probability for an electron  
to tunnel through the barrier even if its energy is less than the potential  
of the barrier.

Two metals must be separated by a space of not more than 10nm.

Tunneling current decreases exponentially with separation between tip and sample.

Experimental methods

the 1  
sample

you want to study 0

a sharp metal 1  
tip

mounted on a 0  
piezoelectric crystal tube to be e  
placed in very close proximity to 0

the sample

a mechanism to control the location  
of the tip in the x

y plane parallel to  
the sample surface

90

feedback

loop to control the  
height of the tip above the sample  
(the z

axis)

Basic Set

up

How to operate?

Raster

the

tip

across

the

surface

at

a

distance

of

0

.

31  
—

7

nm

.

While

the

current

between

the

tip

and

sample

surface

is

measured

And

using

the

current

as

a

feedback

signal

The

tip

moves

over

the

sample

with

its

height

being

adjusted

continuously

to

keep

tunneling

current

constant  
..

Tip

position

is

monitored

to

map

the

surface

topography

of

the

sample

Tip

tungsten or PtIr alloy is used."



## Tip movement in 3

D is controlled by piezoelectric arrays.

The distance between tip and sample is between 0.2 and 0.6 nm which produces a tunneling current of about 0.1

10 nA.

Spatial resolution is 0.01 nm in x and y directions. And about 0.002 nm in z

direction giving a true atomic resolution image & in 3

D.

## Modes of Operation

### Constant current

Constant current is maintained between tip and sample..

vertical position of the tip is changed to maintain a constant current between the two.

### Advantages

Provides excellent surface topographic contrast of the surface atomic contours.

### Disadvantage

Slower scan rate.

### Constant height mode

Distance between tip and sample is kept constant.

Constant tip position results in variations in tunneling current.

### Advantage

faster scan rate possible

Imaging the structure of electrode surface

∞

STM images of the  $\text{Au}(111)$  electrode surface

Concluding remarks

Disadvantage of STM:

1.

Making atomically sharp tips remains something of a dark art.

2.

External and internal vibrations from fans, pumps, machinery, building movements, etc. are big problems.

3.

UHV

STM is not easy to build and to handle

4.

The STM can only scan conductive surfaces or thin nonconductive films and small objects deposited on conductive substrates. It does not work with nonconductive materials, such as glass, rock, etc.

5.

The spatial resolution of STM is fantastic, but the time per scan is typically on the order of seconds, which prevents STM from imaging fast kinetics of electrochemical processes.

STM is one of the most powerful imaging tools with an unprecedented precision.

